# Package 'VisiumIO'

April 28, 2025

Title Import Visium data from the 10X Space Ranger pipeline

Version 1.4.0

Description The package allows users to readily import spatial data obtained from either the 10X website or from the Space Ranger pipeline. Supported formats include tar.gz, h5, and mtx files. Multiple files can be imported at once with \*List type of functions. The package represents data mainly as SpatialExperiment objects.

License Artistic-2.0

**Depends** R (>= 4.4.0), TENxIO

Imports BiocBaseUtils, BiocGenerics, BiocIO (>= 1.15.1), jsonlite, methods, S4Vectors, SingleCellExperiment, SpatialExperiment, SummarizedExperiment

**Suggests** arrow, BiocStyle, data.table, knitr, readr, rmarkdown, tinytest

biocViews Software, Infrastructure, DataImport, SingleCell, Spatial

VignetteBuilder knitr

**Encoding UTF-8** 

**Roxygen** list(markdown = TRUE)

RoxygenNote 7.3.2

BugReports https://github.com/waldronlab/VisiumIO/issues

URL https://github.com/waldronlab/VisiumIO

Collate 'TENxSpatialCSV.R' 'TENxSpatialList-class.R' 'TENxSpatialParquet.R' 'TENxVisium-class.R' 'TENxVisiumList-class.R' 'TENxVisiumHD-class.R' 'VisiumIO-package.R' 'utilities.R'

Date 2025-04-04

git\_url https://git.bioconductor.org/packages/VisiumIO

git\_branch RELEASE\_3\_21

git\_last\_commit fe979fe

git\_last\_commit\_date 2025-04-15

2 VisiumIO-package

Repository Bioconductor 3.21
Date/Publication 2025-04-27
Author Marcel Ramos [aut, cre] (ORCID:
<https: 0000-0002-3242-0582="" orcid.org="">),</https:>
Dario Righelli [aut, ctb],
Helena Crowell [aut, ctb]
Maintainer Marcel Ramos <pre><marcel.ramos@sph.cuny.edu></marcel.ramos@sph.cuny.edu></pre>

# **Contents**

VisiumIO-package																
compareBarcodes																
TENxSpatialCSV-class																
TENxSpatialList-class																
TENxSpatialParquet-cla	SS															
TENxVisium-class								 				 				
TENxVisiumHD-class																
TENxVisiumList-class								 				 				

Index 18

VisiumIO-package

VisiumIO: Import Visium data from the 10X Space Ranger pipeline

# Description

The package allows users to readily import spatial data obtained from either the 10X website or from the Space Ranger pipeline. Supported formats include tar.gz, h5, and mtx files. Multiple files can be imported at once with \*List type of functions. The package represents data mainly as SpatialExperiment objects.

# Author(s)

Maintainer: Marcel Ramos <marcel.ramos@sph.cuny.edu> (ORCID)

Authors:

- Dario Righelli [contributor]
- Helena Crowell [contributor]

#### See Also

Useful links:

- https://github.com/waldronlab/VisiumIO
- Report bugs at https://github.com/waldronlab/VisiumIO/issues

compareBarcodes 3

compareBarcodes

Compare barcodes between raw and filtered data

#### **Description**

This function compares the barcodes between raw and filtered data **depending** on the order of processing. Typically, the "raw" barcodes are compared to the "filtered" ones. The presence of raw barcodes in the filtered data are marked as TRUE in the resulting data.frame.

#### Usage

```
compareBarcodes(
  from_resource,
  to_resource,
  spacerangerOut,
  format = c("mtx", "h5"),
  processing = c("raw", "filtered"),
  ...
)
```

#### **Arguments**

from\_resource character(1) The path to the resource file whose barcodes are used as the basis

of the comparison; typically, the "raw" feature barcodes are used.

to\_resource character(1) The path to the resource file whose barcodes are compared to the

from\_resource; typically, the "filtered" feature barcodes.

spacerangerOut character(1) A single string specifying the path to the sample directory of

spaceranger count. The directory must contain the (processing)\_feature\_bc\_matrix

and spatial sub directories in addition to the outs folder.

format The format of the output. If missing and con is a file name, the format is derived

from the file extension. This argument is unnecessary when con is a derivative

of BiocFile.

processing character(2) A vector of length 2 that corresponds to the processing type. The

processing types are typically "raw" and "filtered". These are the prefixes of the folder names raw\_feature\_bc\_matrix and filtered\_feature\_bc\_matrix. The order of the vector determines the comparison. By default, processing = c("raw", "filtered"), which means barcodes in the raw data are compared

to the filtered data.

... Additional arguments passed to TENxH5 or TENxFileList.

#### Value

A data frame with barcodes of the first element in the processing data type as the first column and a logical vector indicating whether the barcodes are found in the second element in processing. For example, if processing is c("raw", "filtered"), then the first column will be the barcodes in the raw data and the second column will be a logical vector indicating whether the barcodes are found in the filtered data.

#### **Examples**

```
if (interactive()) {
    compareBarcodes(
        from_resource = "V1_Adult_Mouse_Brain_raw_feature_bc_matrix.tar.gz",
        to resource =
            "V1_Adult_Mouse_Brain_filtered_feature_bc_matrix.tar.gz",
   )
   compareBarcodes(
        from_resource =
            "V1_Adult_Mouse_Brain_raw_feature_bc_matrix.h5",
        to_resource =
            "V1_Adult_Mouse_Brain_filtered_feature_bc_matrix.h5"
    )
   compareBarcodes(spacerangerOut = "~/data/outs", format = "h5")
   compareBarcodes(
        spacerangerOut = "~/data/feature_bc_matrix", format = "mtx"
   )
   compareBarcodes(
        spacerangerOut = "~/data/folder_feature_bc_matrix", format = "mtx"
}
```

TENxSpatialCSV-class Represent and import spatial CSV data from 10X Genomics

#### **Description**

TENxSpatialCSV is a class to represent and import spatial CSV files with specific column names. It is a composed class of TENxIO::TENxFile and contains additional slots for the column names and whether the CSV is a list-type of file.

#### Usage

```
TENxSpatialCSV(resource, colnames = .TISSUE_POS_COLS)
## S4 method for signature 'TENxSpatialCSV,ANY,ANY'
import(con, format, text, ...)
```

# Arguments

```
resource character(1) The path to the file

colnames character() A vector specifying the column names of the CSV, defaults to c("barcode", "in_tissue", "array_row", "array_col", "pxl_row_in_fullres", "pxl_col_in_fullres"). Mainly used for the "positions" CSV type of file which does not include column names in the file.
```

con	The connection from which data is loaded or to which data is saved. If this is a character vector, it is assumed to be a file name and a corresponding file connection is created and then closed after exporting the object. If it is a BiocFile derivative, the data is loaded from or saved to the underlying resource. If missing, the function will return the output as a character vector, rather than writing to a connection.
format	The format of the output. If missing and con is a file name, the format is derived from the file extension. This argument is unnecessary when con is a derivative of BiocFile.
text	If con is missing, this can be a character vector directly providing the string data to import.
	Additional inputs to the low level class generator functions

#### **Details**

Typically, the user will not create an object of this class directly but rather use the TENxVisium() constructor function to create an object of this class in the background. The column names are set to the default values of c("barcode", "in\_tissue", "array\_row", "array\_col", "pxl\_row\_in\_fullres", "pxl\_col\_in\_fullres"). The column names can be changed by specifying the colnames argument in the constructor function.

Set the option "VisiumIO.csvreader" to either "data.table" or "readr" to use the data.table::fread or readr::read\_csv functions, respectively. These options are useful when the CSV file is relatively large and the user wants to use faster read-in options. Note that the outputs will still be converted to DataFrame when incorporated to the SpatialExperiment or SingleCellExperiment object.

#### Value

```
TENxSpatialCSV: An object of class TENxSpatialCSV import-method: A DataFrame object containing the data from the CSV file
```

#### **Slots**

```
isList logical(1) A scalar specifying whether the CSV is a list-type of file
colnames character() A vector specifying the column names of the CSV
variant character(1) A scalar specifying the variant of the CSV file "positions", "cell_boundaries",
    or "other". The variant is determined by the name of the CSV file within the constructor func-
tion. Values include "positions", "cell_boundaries", and "other".
compressed logical(1) A scalar specifying whether the CSV is compressed (mainly with a .gz
file extension).
```

```
sample_dir <- system.file(
    file.path("extdata", "10xVisium", "section1"),
    package = "SpatialExperiment"
)</pre>
```

```
spatial_dir <- Filter(
  function(x) endsWith(x, "spatial"), list.dirs(sample_dir)
)
csvresource <- file.path(spatial_dir, "tissue_positions_list.csv")
TENxSpatialCSV(csvresource)
head(import(TENxSpatialCSV(csvresource)), 4)
import(TENxSpatialCSV(csvresource)) |>
  attr("metadata") |>
  lapply(names)
```

TENxSpatialList-class A class to represent and import spatial Visium data

#### **Description**

This class is a composed class of TENxFileList, which can contain a list of TENxFile objects, and a TENxSpatialList object. It is meant to handle spatial Visium data from 10X Genomics.

#### Usage

```
TENxSpatialList(
  resources,
  sample_id = "sample01",
  images = c("lowres", "hires", "detected", "aligned", "aligned_fiducials", "cytassist"),
  jsonFile = .SCALE_JSON_FILE,
  tissuePattern = "tissue_positions.*",
  bin_size = character(0L),
  ...
)

## S4 method for signature 'TENxSpatialList,ANY,ANY'
import(con, format, text, ...)
```

# Arguments

resources	A TENxFileList object or a file path to the tarball containing the matrix / assay data resources.
sample_id	character(1) A single string specifying the sample ID.
images	character() A vector specifying the images to be imported; can be one or multiple of "lowres", "hires", "detected", "aligned".
jsonFile	character(1) A single string specifying the name of the JSON file containing the scale factors.
tissuePattern	character(1) A single string specifying the pattern to match the tissue positions file.
bin_size	character(1) The bin size of the images to import. The default is 008. It corresponds to the directory name square_000um where 000 is the bin value.

... Parameters to pass to the format-specific method.

con The connection from which data is loaded or to which data is saved. If this

is a character vector, it is assumed to be a file name and a corresponding file connection is created and then closed after exporting the object. If it is a BiocFile derivative, the data is loaded from or saved to the underlying resource. If missing, the function will return the output as a character vector, rather than

writing to a connection.

format The format of the output. If missing and con is a file name, the format is derived

from the file extension. This argument is unnecessary when con is a derivative

of BiocFile.

text If con is missing, this can be a character vector directly providing the string data

to import.

#### **Details**

Typically, the user will not create an object of this class directly but rather use the TENxVisium() constructor function to create an object of this class.

#### Value

A SpatialExperiment object

# Methods (by generic)

 import(con = TENxSpatialList, format = ANY, text = ANY): Import a TENxSpatialList object

# Slots

images character() The image name(s) to use with grep and include in the list of files. Can be one of "lowres", "hires", "lowres", "hires", "detected", "aligned", "aligned\_fiducials", or "cytassist".

scaleJSON character(1) The file name of the scale factors JSON file, defaults to 'scalefactors\_json.json'.

tissuePos character(1) The file name of the tissue positions file; typically a .parquet or .csv file

sampleId character(1) A scalar specifying the sample identifier.

binSize The bin size of the images to import. The default slot value is character(). It typically corresponds to the directory name square\_000um where 000 is the bin value.

```
spatial_dir <- system.file(
    file.path("extdata", "10xVisium", "section1", "outs", "spatial"),
    package = "SpatialExperiment"
)

TENxSpatialList(resources = spatial_dir, images = "lowres")</pre>
```

```
TENxSpatialList(resources = spatial_dir, images = "lowres") |>
    metadata() |> lapply(names)
```

TENxSpatialParquet-class

Represent and import spatial Parquet data from 10X Genomics

# **Description**

TENxSpatialParquet is a class to represent and import spatial Parquet files with specific column names. It is a composed class of TENxIO::TENxFile and contains additional slots for the column names and whether the Parquet is a list-type of file.

#### Usage

```
TENxSpatialParquet(resource, colnames = .TISSUE_POS_COLS)
## S4 method for signature 'TENxSpatialParquet,ANY,ANY'
import(con, format, text, ...)
```

#### **Arguments**

resource	character(1) The path to the file
colnames	character() A vector specifying the column names of the Parquet, defaults to c("barcode", "in_tissue", "array_row", "array_col", "pxl_row_in_fullres" pxl_col_in_fullres").
con	The connection from which data is loaded or to which data is saved. If this is a character vector, it is assumed to be a file name and a corresponding file connection is created and then closed after exporting the object. If it is a BiocFile derivative, the data is loaded from or saved to the underlying resource. If missing, the function will return the output as a character vector, rather than writing to a connection.
format	The format of the output. If missing and con is a file name, the format is derived from the file extension. This argument is unnecessary when con is a derivative of BiocFile.
text	If con is missing, this can be a character vector directly providing the string data to import.
	Additional inputs to the low level class generator functions

#### **Details**

Typically, the user will not create an object of this class directly but rather use the TENxVisium() constructor function to create an object of this class in the background. The column names are set to the default values of c("barcode", "in\_tissue", "array\_row", "array\_col", "pxl\_row\_in\_fullres", "pxl\_col\_in\_fullres"). The column names can be changed by specifying the colnames argument in the constructor function.

TENxVisium-class 9

#### Value

```
TENxSpatialParquet(): An object of class TENxSpatialParquet import-method: A DataFrame object containing the data from the Parquet file
```

#### **Slots**

colnames character() A vector specifying the column names of the Parquet

#### **Examples**

```
sample_dir <- system.file(
    file.path("extdata", "binned_outputs", "square_002um", "spatial"),
    package = "VisiumIO"
)
spatial_dir <- Filter(
    function(x) endsWith(x, "spatial"), list.dirs(sample_dir)
)
parquetres <- file.path(spatial_dir, "tissue_positions.parquet")
TENxSpatialParquet(parquetres)
import(TENxSpatialParquet(parquetres))
## metadata in attributes
import(TENxSpatialParquet(parquetres)) |>
    attr("metadata") |>
    lapply(names)
```

TENxVisium-class

A class to represent and import a single Visium Sample

# Description

This class is a composed class of TENxFileList which can contain a list of TENxFile objects and a TENxSpatialList object. It is meant to handle a single Visium sample from 10X Genomics.

# Usage

```
TENxVisium(
  resources,
  spatialResource,
  spacerangerOut,
  sample_id = "sample01",
  processing = c("filtered", "raw"),
  format = c("mtx", "h5"),
  images = c("lowres", "hires", "detected", "aligned"),
  jsonFile = .SCALE_JSON_FILE,
  tissuePattern = "tissue_positions.*\\.csv",
  spatialCoordsNames = c("pxl_col_in_fullres", "pxl_row_in_fullres"),
  ...
```

10 TENxVisium-class

```
## S4 method for signature 'TENxVisium, ANY, ANY'
import(con, format, text, ...)
```

#### **Arguments**

resources A TENxFileList object or a file path to the tarball containing the matrix / assay

data resources.

spatialResource

A TENxSpatialList object or a file path to the tarball containing the spatial data.

spacerangerOut character(1) A single string specifying the path to the sample directory of

spaceranger count. The directory must contain the (processing)\_feature\_bc\_matrix

and spatial sub directories in addition to the outs folder.

sample\_id character(1) A single string specifying the sample ID.

processing character(1) A single string indicating the processing folder available e.g.,

"filtered\_feature\_barcode\_matrix" in the spacerangerOut folder. It can be either "filtered" or "raw" (default "filtered"). Only used when spacerangerOut is

specified.

format The format of the output. If missing and con is a file name, the format is derived

from the file extension. This argument is unnecessary when con is a derivative

of BiocFile.

images character() A vector specifying the images to be imported; can be one or

multiple of "lowres", "hires", "detected", "aligned".

jsonFile character(1) A single string specifying the name of the JSON file containing

the scale factors.

tissuePattern character(1) A single string specifying the pattern to match the tissue posi-

tions file.

spatialCoordsNames

character() A vector of strings specifying the names of the columns in the

spatial data containing the spatial coordinates.

.. In the constructor, additional arguments passed to TENxFileList; otherwise, not

used.

con The connection from which data is loaded or to which data is saved. If this

is a character vector, it is assumed to be a file name and a corresponding file connection is created and then closed after exporting the object. If it is a BiocFile derivative, the data is loaded from or saved to the underlying resource. If missing, the function will return the output as a character vector, rather than

writing to a connection.

text If con is missing, this can be a character vector directly providing the string data

to import.

# Details

Typically, the user will not create an object of this class directly but rather use TENxVisiumList constructor function for multiple samples. Note that the images, jsonFile, tissuePattern, and

TENxVisium-class 11

spatialCoordsNames arguments are only considered when the spacerangerOut argument or both the resources and spatialResource arguments are paths to files.

#### Value

A SpatialExperiment object

#### **Functions**

• import(con = TENxVisium, format = ANY, text = ANY): Import Visium data

#### **Slots**

```
resources A TENxFileList or TENxH5 object containing the Visium data.

spatialList A TENxSpatialList object containing the spatial

coordNames character() A vector specifying the names of the columns in the spatial data containing the spatial coordinates.

sampleId character(1) A scalar specifying the sample identifier.
```

#### See Also

https://support.10xgenomics.com/spatial-gene-expression/software/pipelines/latest/output/overview

```
sample_dir <- system.file(</pre>
    file.path("extdata", "10xVisium", "section1"),
   package = "VisiumIO"
)
## using spacerangerOut folder
tv <- TENxVisium(</pre>
    spacerangerOut = sample_dir, processing = "raw", images = "lowres"
)
import(tv)
## with TENxFileList spacerangerOut input
tvfl <- TENxVisium(</pre>
    spacerangerOut = TENxFileList(sample_dir),
    format = "mtx",
   processing = "raw",
    images = "lowres"
)
import(tvfl)
## check metadata of the object
import(tvfl) |>
   metadata() |>
```

12 TENxVisiumHD-class

```
lapply(names)
## importing h5 format
tvfl <- TENxVisium(</pre>
    spacerangerOut = sample_dir,
    format = "h5",
    processing = "raw",
    images = "lowres"
)
import(tvfl)
rffolder <- file.path(sample_dir, "outs", "raw_feature_bc_matrix")
## using resources and spatialResource inputs
tvfl <- TENxVisium(
    resources = rffolder,
    spatialResource = file.path(dirname(rffolder), "spatial"),
    format = "mtx",
   processing = "raw",
    images = "lowres"
)
import(tvfl)
```

TENxVisiumHD-class

A class to represent and import multiple Visium HD samples

#### **Description**

This class contains a SimpleList of TENxVisiumHD objects each corresponding to one sample. The provided spacerangerOut folder should contain a binned\_outputs folder where multiple bin\_size subfolders are present, e.g., square\_002um.

# Usage

```
TENxVisiumHD(
   resources,
   spatialResource,
   spacerangerOut,
   sample_id = "sample01",
   processing = c("filtered", "raw"),
   format = c("mtx", "h5"),
   images = c("lowres", "hires", "detected", "aligned_fiducials"),
   bin_size = c("008", "016", "002"),
   jsonFile = .SCALE_JSON_FILE,
   tissuePattern = "tissue_positions\\.parquet",
   spatialCoordsNames = c("pxl_col_in_fullres", "pxl_row_in_fullres"),
```

TENxVisiumHD-class 13

```
## S4 method for signature 'TENxVisiumHD, ANY, ANY'
import(con, format, text, ...)
```

#### **Arguments**

resources A TENxFileList object or a file path to the tarball containing the matrix / assay

data resources.

spatialResource

A TENxSpatialList object or a file path to the tarball containing the spatial data.

spacerangerOut character(1) A single string specifying the path to the sample directory of

spaceranger count. The directory must contain the (processing)\_feature\_bc\_matrix

and spatial sub directories in addition to the outs folder.

sample\_id character(1) A single string specifying the sample ID.

processing character(1) A single string indicating the processing folder available e.g.,

"filtered\_feature\_barcode\_matrix" in the spacerangerOut folder. It can be either "filtered" or "raw" (default "filtered"). Only used when spacerangerOut is

specified.

format The format of the output. If missing and con is a file name, the format is derived

from the file extension. This argument is unnecessary when con is a derivative

of BiocFile.

images character() A vector specifying the images to be imported; can be one or

multiple of "lowres", "hires", "detected", "aligned".

bin\_size character(1) The bin size of the images to import. The default is 008. It

corresponds to the directory name square\_000um where 000 is the bin value.

jsonFile character(1) A single string specifying the name of the JSON file containing

the scale factors.

tissuePattern character(1) A single string specifying the pattern to match the tissue posi-

tions file.

spatialCoordsNames

character() A vector of strings specifying the names of the columns in the

spatial data containing the spatial coordinates.

.. In the constructor, additional arguments passed to TENxFileList; otherwise, not

used.

con The connection from which data is loaded or to which data is saved. If this

is a character vector, it is assumed to be a file name and a corresponding file connection is created and then closed after exporting the object. If it is a BiocFile derivative, the data is loaded from or saved to the underlying resource. If missing, the function will return the output as a character vector, rather than

writing to a connection.

text If con is missing, this can be a character vector directly providing the string data

to import.

14 TENxVisiumHD-class

#### **Details**

Typically, the user will provide a path to a directory containing the output of the spaceranger count command. The spaceranger count command outputs a folder containing the "raw" or "filtered" ()\_feature\_bc\_matrix.

#### Value

A SpatialExperiment object

#### **Functions**

• import(con = TENxVisiumHD, format = ANY, text = ANY): Import Visium HD data from multiple bin sizes

```
vdir <- system.file(</pre>
    "extdata", package = "VisiumIO", mustWork = TRUE
## with spacerangerOut folder
TENxVisiumHD(spacerangerOut = vdir, bin_size = "002", images = "lowres")
TENxVisiumHD(spacerangerOut = vdir, bin_size = "002", images = "lowres") |>
    import()
## indicate h5 format
TENxVisiumHD(
    spacerangerOut = vdir, bin_size = "002",
    images = "lowres", format = "h5"
)
TENxVisiumHD(
    spacerangerOut = vdir, bin_size = "002",
    images = "lowres", format = "h5"
) |>
    import()
## use resources and spatialResource arguments as file paths
TENxVisiumHD(
    resources = file.path(
        vdir, "binned_outputs", "square_002um",
        "filtered_feature_bc_matrix.h5"
   ),
    spatialResource = file.path(
        vdir, "binned_outputs", "square_002um",
        "spatial"
   ),
    bin_size = "002", processing = "filtered",
    images = "lowres", format = "h5"
) |>
    import()
```

TENxVisiumList-class 15

```
## provide the spatialResource argument as a TENxFileList
TENxVisiumHD(
    resources = file.path(
        vdir, "binned_outputs", "square_002um",
        "filtered_feature_bc_matrix.h5"
    ),
    spatialResource = TENxFileList(
        file.path(
            vdir, "binned_outputs", "square_002um",
            "spatial"
        )
     ),
     bin_size = "002", images = "lowres", format = "h5"
) |>
    import()
```

TENxVisiumList-class A class to represent and import multiple Visium samples

# **Description**

This class contains a SimpleList of TENxVisium objects each corresponding to one sample.

# Usage

```
TENxVisiumList(
  sampleFolders,
  sample_ids,
  processing = c("filtered", "raw"),
  images = c("lowres", "hires", "detected", "aligned"),
  format = c("mtx", "h5"),
  jsonFile = .SCALE_JSON_FILE,
  tissuePattern = "tissue_positions.*\\.csv",
  spatialCoordsNames = c("pxl_col_in_fullres", "pxl_row_in_fullres"),
  ...
)

## S4 method for signature 'TENxVisiumList,ANY,ANY'
import(con, format, text, ...)
```

#### Arguments

sampleFolders character() A vector of strings specifying the directories containing the output of the spaceranger count command.

sample\_ids character() A vector of strings specifying the sample IDs. If not provided, the sample IDs will be the names of the sampleFolders. Therefore, the sample\_ids must be the same length as sampleFolders.

16 TENxVisiumList-class

processing character(1) A single string indicating the processing folder available e.g.,

"filtered\_feature\_barcode\_matrix" in the spacerangerOut folder. It can be either "filtered" or "raw" (default "filtered"). Only used when spacerangerOut is

specified.

images character() A vector specifying the images to be imported; can be one or

multiple of "lowres", "hires", "detected", "aligned".

format The format of the output. If missing and con is a file name, the format is derived

from the file extension. This argument is unnecessary when con is a derivative

of BiocFile.

jsonFile character(1) A single string specifying the name of the JSON file containing

the scale factors.

tissuePattern character(1) A single string specifying the pattern to match the tissue posi-

tions file.

spatialCoordsNames

character() A vector of strings specifying the names of the columns in the

spatial data containing the spatial coordinates.

... In the constructor, additional arguments passed to TENxFileList; otherwise, not

used.

con The connection from which data is loaded or to which data is saved. If this

is a character vector, it is assumed to be a file name and a corresponding file connection is created and then closed after exporting the object. If it is a BiocFile derivative, the data is loaded from or saved to the underlying resource. If missing, the function will return the output as a character vector, rather than

writing to a connection.

text If con is missing, this can be a character vector directly providing the string data

to import.

# Details

Typically, the user will provide a path to a directory containing the output of the spaceranger count command. The spaceranger count command outputs a folder containing the "raw" or "filtered" ()\_feature\_bc\_matrix.

# Value

A SpatialExperiment object

#### **Functions**

• import(con = TENxVisiumList, format = ANY, text = ANY): Import multiple Visium samples

#### See Also

https://support.10xgenomics.com/spatial-gene-expression/software/pipelines/latest/output/overview

TENx Visium List-class 17

```
sample_dirs <- list.dirs(
    system.file(
        file.path("extdata", "10xVisium"),
        package = "VisiumIO"
    ),
    recursive = FALSE, full.names = TRUE
)

tvl <- TENxVisiumList(
    sampleFolders = sample_dirs,
    sample_ids = c("sample01", "sample02"),
    processing = "raw",
    images = "lowres",
    format = "mtx"
)

import(tvl)</pre>
```

# **Index**

* internal	TENxSpatialList
VisiumIO-package, 2	(TENxSpatialList-class), 6
.TENxSpatialCSV (TENxSpatialCSV-class),	TENxSpatialList-class, 6
4	TENxSpatialParquet, 9
.TENxSpatialList	TENxSpatialParquet
(TENxSpatialList-class), 6	(TENxSpatialParquet-class), 8
.TENxSpatialParquet	TENxSpatialParquet-class, 8
(TENxSpatialParquet-class), 8	TENxVisium, 15
.TENxVisium (TENxVisium-class), 9	TENxVisium (TENxVisium-class), 9
.TENxVisiumHD (TENxVisiumHD-class), 12	TENxVisium(), 5, 7, 8
.TENxVisiumList (TENxVisiumList-class),	TENxVisium-class, 9
15	TENxVisiumHD, <i>12</i>
	TENxVisiumHD (TENxVisiumHD-class), 12
BiocFile, 3, 5, 7, 8, 10, 13, 16	TENxVisiumHD-class, 12
	TENxVisiumList, <i>10</i>
compareBarcodes, 3	TENxVisiumList (TENxVisiumList-class)
	15
<pre>import,TENxSpatialCSV,ANY,ANY-method</pre>	TENxVisiumList-class, 15
(TENxSpatialCSV-class), 4	
<pre>import,TENxSpatialList,ANY,ANY-method</pre>	VisiumIO(VisiumIO-package), 2
(TENxSpatialList-class), 6	VisiumIO-package, 2
<pre>import,TENxSpatialParquet,ANY,ANY-method</pre>	
(TENxSpatialParquet-class), $8$	
<pre>import,TENxVisium,ANY,ANY-method</pre>	
(TENxVisium-class), 9	
<pre>import,TENxVisiumHD,ANY,ANY-method</pre>	
(TENxVisiumHD-class), 12	
<pre>import,TENxVisiumList,ANY,ANY-method</pre>	
(TENxVisiumList-class), 15	
SpatialExperiment, 11, 14, 16	
TENxFile, 6, 9	
TENxFileList, 6, 9-11, 13, 16	
TENxH5, 11	
TENxIO::TENxFile, 4, 8	
TENxSpatialCSV, 5	
TENxSpatialCSV (TENxSpatialCSV-class), 4	
TENxSpatialCSV-class, 4	
TENxSpatialList. 6, 9-11, 13	