canine2probe

April 29, 2025

| canine2probe Probe sequence for microarrays of type canine2. | canine2probe | Probe sequence for microarrays of type canine2. | |
|--|--------------|---|--|
|--|--------------|---|--|

Description

This data object was automatically created by the package AnnotationForge version 1.11.21.

Usage

```
data(canine2probe)
```

Format

A data frame with 473162 rows and 6 columns, as follows.

| sequence | character | probe sequence |
|------------------------------|-----------|-------------------------------------|
| X | integer | x-coordinate on the array |
| y | integer | y-coordinate on the array |
| Probe.Set.Name | character | Affymetrix Probe Set Name |
| Probe.Interrogation.Position | integer | Probe Interrogation Position |
| Target.Strandedness | factor | Target Strandedness |

Source

The probe sequence data was obtained from http://www.affymetrix.com. The file name was Canine $\2\$ probe $\$ tab.

Examples

```
canine2probe
as.data.frame(canine2probe[1:3,])
```

Index

* datasets
canine2probe, 1
canine2probe, 1